

aidPlus

Orthopedic Supports



A core step towards health!

A brand presented by **AbdulWahed Bin Shabib Investment Group LLC.**



Our orthopedic support and braces are made specifically to address musculature problems and are used to support, align, stabilize, and protect specific body components, including the muscles, joints, and bones as they recover from trauma or injury.



INDEX

1.	Neck 1.1 1.2 1.3 1.4	Collar Soft Collar Philadelphia Collar Cervical Collar with Chin Emergency Collar	01
2.	Arm-\$ 2.1 2.2 2.3 2.4	Shoulder Support Shoulder Support Arm Sling (Luz) Arm Sling (Mesh) Arm Sling (Simple)	03
3.	Elbow 3.1 3.2 3.3	Support Elastical Elbow Support Padded Elbow Support Knitted Epicondylitis Elbow Support	05
4.	Hand 4.1 4.2 4.3 4.4 4.5	Splint Hand Wrist Splint (Both Sides) Hand & Wrist Splint (Right-Left) Hand Wrist Splint with Thumb Grip (Rig Neoprene Wrist Splint (Right-Left) Neoprene Wrist Splint with Thumb Grip (Right-Left)	07 ht-Left)
5.	Finge 5.1 5.2 5.3	r Splint Flat Finger Baseball Splint Frog Splint	09
6.	Corse 6.1 6.2 6.3 6.4 6.5 6.6	wit Support Woolen Corset Woolen Corset (with Velcro) Camel Fur Corset Dog Fur Corset Abdominal Corset Lumbosacral Corset with Belt (26 cm/3) Maternity Corset Belt	11 2 cm, Black/White)
7.	Hernia	a Ties (Double)	13
8.	Surgio 8.1 8.2 8.3	cal Stocking CCL-1 Under Knee Surgical Stocking CCL-1 Pregnancy Surgical Stocking CCL-2 Over Knee Surgical Stocking	14

9.	Knee Support 9.1 Wool Knee Brace 9.2 Closed Knee Support 9.3 Patella Knee Support 9.4 Patella & Ligament Knee Support 9.5 Knee Support (Steel Hinged) 9.6 Closed Patella Knee Support 9.7 Immobilizer	15
10.	Knee Braces 10.1 Immobilizer	17
11.	Thigh & Calf Support 10.1 Thigh Support 10.2 Calf Support	18
11.	Ankle Support 11.1 Eight Bandage 11.2 Cross Strap Ankle Support 11.3 Ankle with Ligament Support	19
12.	Foot Health Product 12.1 Silicone Insoles 12.2 Silicone Insoles (Ark High) 12.3 Silicone Insoles (with Textile) 12.4 Silicone Heel-cups 12.5 Silicone Yoyo Heel-cups 12.6 Silicone Toe Separator 12.6 Silicone Metatarsal Pad 12.7 Silicone Cast Shoes	20
13.	Life Support 13.1 Thermojel 13.2 Hemorrhoid Bathtub	22
14.	One Size Product 14.1 Dequrvain Splint 14.2 Bandage Wrist Support 14.3 Simple Wrist Support 14.4 Patella Tendon Band 14.5 Lumbosacral Corset	23
15.	Elastic Support 15.1 Elastic Elbow Support 15.2 Elastic Knee Support 15.3 Elastic Ankle Support	24
16.	Visco & Polyurethane 16.1 Visco Snoring Pillow 16.2 Visco Neck Pillow 16.3 Lumbar Cushion 16.4 Seat Cushion (Open) 16.5 Travel Pillow	25

Neck Collar

We offer a range of neck collars, depending upon the type of neck injury you are experiencing or the underlying cause of your neck pain.

Neck collars are a common form of treatment for some neck injuries, surgeries, and pain. Its applications include cervical hernia, strained muscles, osteoarthritis, and muscular spasms.



Soft Collar (Sponge)

It is designed according to the anatomical neck shape, the sponge is a soft cervical collar. For patients, to avoid excessive neck motion, it allows recalling sensory feedback and keeps warm around the neck.

Sizes		Q
Small	Medium	Large
28 cm	34 cm	40 cm
35 cm	39 cm	45 cm



Philadelphia Collar

A Philadelphia neck collar is a stiff, latex-free, non-toxic, hypoallergenic, foam collar that supports your neck in the front and behind with rigid plastic. Restricting any movement while healing occurs, is intended to stabilize the cervical spine.

Sizes		Q
Small	Medium	Large
33 cm	38 cm	42 cm
37 cm	42 cm	46 cm



Cervical Collar with Chin

A cervical collar is used to support the neck. Cervical collars are also used therapeutically, to help stabilize the cervical column and relieve pain. It provides chin support by enhancing immobilization with a chin lock. Hard collars are most commonly prescribed for cervical spine stabilization following either trauma, surgery, or fractures/dislocations.

Sizes				3
Small	Medium	Large	XL	XXL
32 cm	36 cm	40 cm	43 cm	46 cm
35 cm	39 cm	43 cm	46 cm	50 cm

Neck Collar





Emergency Collar

1.5 mm, plastic made, eva coated inner surface emergency collar is used as a first cervical collar in emergency situations.

Standard Size





Arm-Shoulder Support



Shoulder Supporter

Made of the highly resilient substance neoprene. Prevents edema by using the compression effect. Reduces muscle tone, which lessens discomfort. Underarm Velcro straps provide comfort and convenience of use.

Application

- Semi-dislocations
- Tendinitis
 Osteoarthritis

Used after surgery

Small	Mediun	ո Large	XL	XXL
22 cm	24 cm	27 cm	30cm	34 cm
34 cm	27 cm	30 cm	34 cm	38 cm

Elbow Support





Arm Sling (Lux)

The inner surface is made of towel, the exterior is made of sponge like cloth and the belt provides extra protection and support.

Application

- Shoulder dislocations
- Postoperative resting
- Preventing fractures of distal edema formation

Sizes

Small	Medium	Large	XL	XXL
30 cm	33 cm	36 cm	39 cm	42 cm
33 cm	36 cm	39 cm	42 cm	45 cm



Arm Sling (Mesh)

Inner and outer surface is made up of breathable perforated fabric which allows the air to pass for preventing the sweat and discomfort.

Application

Shoulder dislocations

- Postoperative resting
- Preventing fractures of distal edema formation

Sizes

Small	Medium	Large	XL	XXL
30 cm	33 cm	36 cm	39 cm	42 cm
33 cm	36 cm	39 cm	42 cm	45 cm



Arm Sling (Simple)

The inner surface is made up of towel like material while the outer surface is made up of sponge like fabric to provide extra comfort.

Size

Adjustable measurement

Elbow Support



Elastical Elbow Support

Made by neoprene material, designed according to anatomy.

Application

- · Protects against trauma
- · For heat and compression
- Protects during and after sports activities.

Small	Medium	Large	XL	XXL	
22 cm	24 cm	26 cm	28cm	30 cm	
24 cm	26 cm	28 cm	30 cm	32 cm	







Padded Elbow Support

Padded Elbow Support is made of neoprene material and is anatomically designed. The elbow part is supported with a specially designed eva material pad, providing more protection.

Application

- · Protects against trauma
- Olecranon bursitis
- Protects during and after sports and physical activities.

Small	Medium	Large	XL	XXL
22 cm	24 cm	26 cm	28cm	30 cm
24 cm	26 cm	28 cm	30 cm	32 cm



Knitted Epicondylitis Elbow Support

Produced using neoprene materials. The anatomically developed silicon pad tendons' droplets that attach to the muscle groups and compress the lateral epicondyle relieve stress and pain.

Application

- · Epicondylitis lateral
- · Increase in the joint fluid
- In cases of elbow pain.

Small	Medium	Large	XL	XXL	
22 cm	24 cm	26 cm	28cm	30 cm	
24 cm	26 cm	28 cm	30 cm	32 cm	

Hand Splint

Used to protect and support the dislocated or broken bones and joints of the hand. Carpal tunnel syndrome, tendonitis, tenosynovitis, wrist pain, and protection following injury and surgery are the principal applications.





Made from breathable cotton fabric. In order to offer immobilization and support, the anatomically formed wrist support is made of aluminum volar. The wrist stabilization is protected by a bandage on it.

Sizes

Small	Medium	Large	XL
15 cm	19 cm	19 cm	23 cm



Hand & Wrist Splint (Right-Left)

Provides immobilization and constrains to the movements of the wrist with the help of aluminum rod inside.

Sizes

Small	Medium	Large	XL	XXL
15 cm	17 cm	19 cm	21 cm	23 cm
17 cm	19 cm	21 cm	23 cm	25 cm

07

Hand Splint





Hand Wrist Splint

with Thumb Grip (Right-Left)

Made of breathable cotton fabric that has anatomically contoured with aluminum support for restricting the movement of wrist and thumb. The bandage around the ankle is provided stability.

-		
C .	70	•
-31	/-	

Small	Medium	Large	XL	XXL	
15 cm	17 cm	19 cm	21 cm	23 cm	
17 cm	19 cm	21 cm	23 cm	25 cm	



Neoprene Wrist Splint (Right-Left)

Neoprene Wrist Splint, made of neoprene material. Anatomically shaped aluminum support provides immobilization, shaping natural wrist position. Increases stabilization around the wrist.

Small	Medium	Large	XL	XXL	
15 cm	17 cm	19 cm	21 cm	23 cm	
17 cm	19 cm	21 cm	23 cm	25 cm	



Neoprene Wrist Splint

with Thumb Grip (Right-Left)

Is made of neoprene material and has got a anatomically contoured aluminum support to restrict the movement of wrist and thumb. The bandage around the wrist provides stability.

Sizes

Small	Medium	Large	XL	XXL
15 cm	17 cm	19 cm	21 cm	23 cm
17 cm	19 cm	21 cm	23 cm	25 cm

Finger Splint

The finger is immobilized by the extension of splinter that has unique aluminum coating on its inside surface. Moreover, apply to finger fractures, sprains, and injuries to the tendon and nerve.



Corset Support







Woolen Corset (with Velcro)

Wool fabric is used in corset support production, providing warmth and compression in the waist area.

Sizes

Small	Medium	Large	XL	XXL	
75 cm	85 cm	95 cm	105 cm	115 cm	
85 cm	95 cm	105 cm	115 cm	125 cm	



Camel Fur Corset

It is constructed with a fabric similar to rabbit fur that is knitted using a unique technique. It can be applied to relieve rheumatic aches as well as to keep the waist, belly, and back warm.

Small	Medium	Large	XL	XXL
75 cm	85 cm	95 cm	105 cm	115 cm
85 cm	95 cm	105 cm	115 cm	125 cm

Corset Support



Dog Fur Corset

It is made with a fabric that resembles dog hair and a special knitting technique. It can be used to keep the back, waist, and belly warm as well as to ease rheumatic pains.

Sizes

Small	Medium	Large	XL	XXL
75 cm	85 cm	95 cm	105 cm	115 cm
85 cm	95 cm	105 cm	115 cm	125 cm



Abdominal Corset

It is a corset that supports the abdomen and is made of flexible, breathable fabric that prevents sweating.

Sizes

Small	Medium	Large	XL	XXL
75 cm	85 cm	95 cm	105 cm	115 cm
85 cm	95 cm	105 cm	115 cm	125 cm





Lumbosacral Corset with Belt

(26 cm/32 cm, Black/White)

Produced from breathable fabric that can reduce the load in disc, thanks to its extra cross wise bands.

Small	Medium	Large	XL	XXL	
75 cm	85 cm	95 cm	105 cm	115 cm	
85 cm	95 cm	105 cm	115 cm	125 cm	

Corset Support





Maternity Corset Belt

Composed of cotton fabric, which is elastic and breathable. Maternity Corset Belt can be worn during pregnancy and after delivery. Secure supports is brought to the front part of the belly in the transition of daily activities.

Small	Medium	Large	XL	XXL
75 cm	85 cm	95 cm	105 cm	115 cm
85 cm	95 cm	105 cm	115 cm	125 cm

Hernia Ties (Double)

Made of supple, stretchable cotton fabric, that intended to push back hernia. It can be utilized before and after surgery to stop recurrence and aid in healing. It provides the ideal support, comfort, and prevention of pain and suffering. It is simple and convenient to utilize in daily life. Can be worn underneath clothing.



Surgical Stocking




After surgery, surgical stockings are used to avoid deep vein thrombosis (DVT), a condition where blood clots form in the leg. Indications include varicosities, lymphedema, venous eczema and ulceration, deep vein thrombosis, and post thrombotic syndrome in the lower extremities due to venous and lymphatic insufficiency.



CCL-1 Under Knee Surgical Stocking



CCL-2 Over Knee Surgical Stocking



CCL-1 Pregnancy Surgical Stocking

Knee Support

Its perfect for knees that are stiff, weak, or sore. Firm support is provided by a flexible and lightweight construction without limiting mobility. Keeps your knee mobile while assisting in the relief of discomfort and swelling. Additionally, provides safety before, during, and after sports or physical activity. The contour design ensures a perfect fit by preventing bunching and slipping.



Small	Medium	Large	XL	XXL
31 cm	34 cm	37 cm	40 cm	44 cm
34 cm	37 cm	40 cm	44 cm	48 cm

Patella Knee Support

Knee Support



Patella & Ligament Knee Support



Closed Patella Knee Support

Knee Support (Steel Hinged)

Small	Medium	Large	XL	XXL
31 cm	34 cm	37 cm	40 cm	44 cm
34 cm	37 cm	40 cm	44 cm	48 cm

Knee Braces



Immobilizer

It is a long, fabric brace that wraps around the leg and shin. On the interior and outside of the brace are metal support struts. The knee immobilizer is secured in place with Velcro straps. The interior has a spongy cloth towel covering it. Use for postoperative meniscal injuries and relaxing post surgery.

Small	Medium	Large
22 cm	24 cm	26 cm
24 cm	26 cm	28 cm

Thigh & Calf Support

Supports the thigh and calf to relieve pain, give heat, compression, and muscular stress before, during, and



Ankle Support

These foot compression sleeves will offer much-needed all-natural pain alleviation for foot ailment, including plantar fasciitis. In addition, they offer ankle and heel arch support for people who have chronic illnesses or injuries that prevent them from fully enjoying life. For added arch support, this soft, light sock can be worn underneath regular socks or athletic socks. They aid blood flow to lessen fatigue and encourage post-exercise recovery.



Eight Bandage

Cross Strap Ankle Support

Sizes

Small	Medium	Large	XL	XXL	
31 cm	34 cm	37 cm	40 cm	44 cm	
34 cm	37 cm	40 cm	44 cm	48 cm	



Ankle Stabilizer with Ligament Support

	Small	Medium	Large	XL	XXL	
	21 cm	23 cm	25 cm	27 cm	29 cm	
9	23 cm	25 cm	27 cm	29 cm	31 cm	



Foot Health Product

Improved foot and leg alignment, increased comfort, and relief from stress and discomfort brought on by flat feet, bunions, arthritis, plantar fasciitis, and Achilles tendonitis are all benefits of arch support. Keeps your feet in the proper position. Stabilizes the foot and shields your heel from the forceful pressure of landing whether you're running or walking.



Foot Health Product



Silicone Heel-cups



Silicone Yoyo Heel-cups

Sizes

Small	Medium	Large
35 cm	38 cm	42 cm
38 cm	41 cm	45 cm



Silicone Toe Separator

Sizes

Small	Medium	Large
34 cm	38 cm	42 cm
37 cm	41 cm	45 cm



Metatarsal Silicone Pad

Standard Size



Cast Shoe

Small	Medium	Large
33 cm	36 cm	41 cm
35 cm	40 cm	45 cm

Life Support



Hemorrhoid Bathtub

One Size Products



Elastic Support



Sizes

Small	Medium	Large	XL
22 cm	24 cm	26 cm	28 cm
24 cm	26 cm	28 cm	30 cm



Sizes

Small	Medium	Large	XL
31 cm	34 cm	37 cm	40 cm
34 cm	37 cm	40 cm	44 cm

Elastic Knee Support



Elastic Ankle Support

Small	Medium	Large	XL	
19 cm	21 cm	23 cm	25 cm	
21 cm	23 cm	25 cm	27 cm	

Visco & Polyurethane



Visco Snoring Pillow

Elastic Visco is used in its production. It aids in the prevention of snoring. It reduces the headache and dry mouth carried on by snoring. It makes breathing, heavy nighttime perspiration, and restless sleep easier.



Visco Neck Pillow

The ergonomic design of the Visco Neck Pillow conforms to the curves of your head and neck to maintain proper alignment throughout the night during sleep, preventing pain. Stiffness and pain are lessened by improved body alignment, particularly in the head, neck, and shoulders.




Visco & Polyurethane




Lumbar Cushion

This lumbar cushion assists you ease on your office chairs while maintaining the natural S curve of your spine. It also supports your core muscles. It eases stress in the upper, middle, and lower bodies. This chair backrest pillow's perfect construction makes it suitable for usage both at home and in the office.



Seat Cushion (Open)

The seat cushion has a soft, pleasant feel because it adjusts to the curve of your body. Provides the optimum support for your lower back. Reduces tension and soreness in the lower back caused by prolonged sitting.



Travel Pillow

This comfy pillow allows restful sleeps while traveling, making it ideal for those long flight and car travel. The distinctive design is curved for superior neck support, and the incredibly soft flocked surface offers comfort.





Contact Us

Phone:

+971 4210 28 92

Mail:

info@bshs.ae

Address:

Shop #17, Bin Shabib Mall, Baghdad Street Al Qusais, Dubai-UAE, Zip- 11370